

## TTC 0.05% SOLUTION

### VIAL CONTENTS

Triphenyltetrazolium chloride 0.0125 g  
Distilled water 25 mL

### DIRECTIONS

The solution is ready to use. Prepare 500 mL of TTC Tergitol 7 Agar Base (REF 402160T) sterilised and cooled to 50°C and add the contents of one vial of TTC 0.05 % Solution (25 mL). Mix well and distribute.

### INTENDED USE

For the isolation and presumptive identification of coliform bacteria by means of the reduction of TTC to formazan, producing pink-red colonies. This supplement is recommended by ISO 7899-2. For the details of the procedure see the technical sheets of TTC Tergitol 7 Agar (REF 402160T)

### STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

### PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

### PACKAGING

**421510 TTC, 0.05% Solution, 10 x 25 mL vials (each vial for 500 mL of medium)**