

YERSINIA SELECTIVE SUPPLEMENT

VIAL CONTENTS (for 500 ml of medium)

Cefsulodin 7.50 mg

Novobiocin 1,25 mg

DIRECTIONS

Aseptically reconstitute one vial with 5 ml of sterile distilled water. Prepare 500 ml of Cefsulodin-Irgasan-Novobiocin (CIN) Agar Base (cat. N° 401302), autoclaved and cooled to 50°C. Add the supplement, mix well and distribute into sterile Petri dishes.

INTENDED USE

For the selective isolation of *Yersinia enterocolitica* from specimens of clinical, food or environmental origin. For the details of the procedure see the technical sheet of CIN Agar Base, cat. N° 401302.

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory. Consult the material safety data sheet before the use.

PACKAGING**4240011****Yersinia Selective Supplement,****10 vials each for 500 ml of medium**